



optimum medical  
support™

## Training and education to Optimise patient care

Visit [www.optimummedicalsupport.co.uk](http://www.optimummedicalsupport.co.uk) to find out more.



Informative user guide available to assist clinicians to meet best practice guidelines on patient and carer education.

Email [ugo@optimummedical.co.uk](mailto:ugo@optimummedical.co.uk) to request your copy.

## REFERENCES

- <sup>1</sup> Quality Governance in the NHS (2011) [https://www.gov.uk/government/uploads/system/uploads/attachment\\_data/file/216321/dh\\_125239.pdf](https://www.gov.uk/government/uploads/system/uploads/attachment_data/file/216321/dh_125239.pdf)
- <sup>2</sup> National Institute for Health and Care Excellence (2012) Guidance 139. Infection: Prevention and Control of Healthcare-associated Infections in Primary and Community Care. London: National Institute for Health and Care Excellence; Available at HYPERLINK "http://www.nice.org.uk" www.nice.org.uk.
- <sup>3</sup> Loveday HP, Wilson JA, Pratt RJ et al. (2013) Epic3: National Evidence-Based guidelines for Preventing Healthcare-Associated Infections in NHS Hospitals in England. Journal of Hospital Infection; 86S1: S1-S70.



optimummedical

### Optimum Medical

First floor, Unit 5, Armley Court, Armley Road,  
Leeds, LS12 2LB, United Kingdom  
Tel.: +44 (0) 8456 435 479, Fax: +44 (0) 1142 383 826  
[www.optimummedical.co.uk](http://www.optimummedical.co.uk)



Company Registration No: 6779883 OMS-UDB-102 Sept 15

independence with confidence

**ugo™**  
leg bags

A THOUGHTFUL RANGE OF URINARY LEG BAGS BASED ON PATIENT USER EVIDENCE

Designed to be used with either an indwelling urethral or suprapubic catheter



optimummedical®

[www.optimummedical.co.uk](http://www.optimummedical.co.uk)



independence  
with confidence



At **Optimum Medical** we believe we have a responsibility to provide enough information that those using our products have a **safe, effective and positive experience\***.

**Ugo Leg Bags** are designed with a range of features to meet individual patient and clinical requirements\*\*, and each is provided with a comprehensive patient and carer user guide to make the experience using our products the best it can be.

\* Darzi's Definition of High Quality Care<sup>1</sup>  
\*\* In line with NICE<sup>2</sup> and EPIC<sup>3</sup> best practice guidelines

- ← All Ugo Leg Bags are sterile pouch packaged which is printed with key user information in a clear and simple format
- All Ugo Leg Bags can remain attached to an indwelling catheter for up to 7 days

### COMMUNITY LEG BAGS

Quick dry fabric backing for patient comfort  
1 pair of **Ugo Fix Leg Bag Straps** in every box

Pr. Name	Pr. Code	Bag Vol.	Tube Length	Tap Design	Inner Box	Outer Case
Ugo 1C	1C/500/ST/L	500ml	Short Tube	Lever Tap	Box of 10	8 Boxes of 10
Ugo 2C	2C/500/LT/L	500ml	Long Tube	Lever Tap	Box of 10	8 Boxes of 10
Ugo 3C	3C/500/ST/T	500ml	Short Tube	T tap	Box of 10	8 Boxes of 10
Ugo 4C	4C/500/LT/T	500ml	Long Tube	T tap	Box of 10	8 Boxes of 10
Ugo 5C	5C/350/DI/L	350ml	Direct Inlet	Lever Tap	Box of 10	8 Boxes of 10
Ugo 6C	6C/350/ST/L	350ml	Short Tube	Lever Tap	Box of 10	8 Boxes of 10
Ugo 7C	7C/750/ST/L	750ml	Short Tube	Lever Tap	Box of 10	8 Boxes of 10
Ugo 8C	8C/750/LT/L	750ml	Long Tube	Lever Tap	Box of 10	8 Boxes of 10

### HOSPITAL LEG BAGS

Non fabric backing to meet infection control standards in clinical environments 1 pair of **Ugo Fix Leg Bag Straps** in every pouch

Ugo 1H	1H/500/ST/L	500ml	Short Tube	Lever Tap	Box of 10	8 Boxes of 10
Ugo 2H	2H/500/LT/L	500ml	Long Tube	Lever Tap	Box of 10	8 Boxes of 10
Ugo 3H	3H/500/ST/T	500ml	Short Tube	T tap	Box of 10	8 Boxes of 10
Ugo 4H	4H/500/LT/T	500ml	Long Tube	T tap	Box of 10	8 Boxes of 10

### LEG BAG STRAPS

Pr. Name	Pr. Code	Inner Box	Outer Case
Ugo Fix Leg Bag Straps	3011	10 Pairs*	20 Boxes of 10

\* Each pair contains one 75cm and one 45cm leg bag strap

### FIXATION DEVICES

To secure and support catheter and drainage systems

The range includes the **Ugo Fix Gentle** catheter clip, **Ugo Fix Catheter Strap** and **Ugo Fix Sleeve**.

Find out more:  
[www.ugofix.co.uk](http://www.ugofix.co.uk)

- Natural leg shape design
- Choice of fabric or non fabric backing depending on clinical environment and patient choice
- Durable but soft medical grade plastic material for strength and reduced noise when moving

LATEX FREE



Comprehensive **Instructions for Use** on how to connect, empty and disconnect bags in a range of European and International languages

## KEY STRENGTHS OF THE RANGE



- 1 Transparent universal stepped connector for secure fitting into all indwelling catheters
- 2 Needle free sprung sample port which locks syringe for safe sample collection
- 3 • Soft kink resistant inlet tubing purpose designed in diameter to reduce bulk and linear pressure to skin  
• Ring welded inlet tubing for connection strength to ensure leak free continuous drainage
- 4 Non return valve to prevent urine returning up the drainage tube
- 5 4 strengthened eyelets to ensure level fixation using enclosed **Ugo Fix Leg Bag Straps**
- 6 Discreet bag printing with unique naming and numbering system for quick and easy re-ordering
- 7 Purpose measured soft silicone outlet tubing for perfect kink-free link drainage connection

Available in:

- 8 Lever tap or 9 T tap

